



# RUDOLF

BECAUSE WE CARE



## NEW

### Piraña Cut LAP Forceps

3 Functions all in one

# Multifunctional Forceps, Piraña Cut



### Scissors

Sharp scissor for the final cut of the suturing process

### Needle holder

Fine structure for a secure grip of the needle during suturing

### Special Coating

The black coating of the insert allows for a longer lifetime of the scissor blades

### Forceps

Teeth 1:2 enable a tight grip of the tissue prior to suturing



5.0  
mm

A complete Screw-In Instrument consists of:  
- Handle, Tube, Insert



Multifunctional Forceps Piraña Cut,  
1:2 teeth, double action, Ø 5 mm

Insert, only	Length	Ordering no.		
RS000-595	330 mm	RS115-595	RS135-595	RS225-595
RS000-795	450 mm	RS117-795	RS127-795	RS227-795